

TK-Snap und TK-Fric System

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TK-Fric

1. Processing techniques

Fabrication of the primary telescope crown


1 For the TK-Snap system, the wax is scraped back on the primary crown in the approximal area, the primary part is attached using the parallel holder and waxed up flush with the surface. The transition from the occlusal to the approximal surface must be rounded in the area of the primary part to enable the friction element to slide over easily.

2 The primary telescope crown can now be invested and cast. The retaining pins of the primary part can remain on the crown, as it is easier to remove them later in metal. This bore is not required for the TK-Fric system, but is advantageous in terms of the compatibility of both systems..

3 The engagement depth of the TK-Fric element is 0.3 mm. If this dimension is undershot due to subsequent milling, the bore depth can be restored using the TKR milling cutter.

Secondary crown in the lift-off technique


4a During the wax-up of the secondary crown, a window is left clear in the area of the snap bore. In this window, the secondary part is positioned exactly over the bore of the primary part using the centring device.

The secondary part is waxed up into the secondary crown and covered from the outside with a wax plate. This prevents liquid wax from flowing into the secondary part in an uncontrolled manner.

5a After casting, all investment material residues must be removed from the secondary part. Only after completion is the TK-Snap or the TK-Fric element inserted into the secondary part. Secondary crown in the one-piece casting technique

Secondary crown in the one-piece casting technique


4b The placeholder is fixed to the primary crown with wax, with the centring button positioned exactly in the bore of the primary crown. For TK-Fric work without a bore in the primary crown, the back of the placeholder lies tangentially against the primary crown.

5b The duplicate model is fabricated from investment material in the usual manner. During this process, the edges in the area of the placeholder must not be damaged.

6b The model casting is waxed up together with the crowns, and the placeholder is integrated into this process.

7b When finishing the model casting, care must be taken to ensure that no investment material residues remain in the receiving cavity for the secondary part. These can prevent the exact fit of the secondary part.

8b Should difficulties in fitting occur when inserting the secondary part into the receiving cavity, the edges of the secondary part may be slightly rounded using a rubber polisher. The bonding of the secondary part can take place either before or after completion. The secondary part is bonded using the proven Si-tec adhesive.

9 Finally, the TK-Snap or TK-Fric friction element is inserted into the secondary part. The friction elements are not bonded, as they are held in the secondary part by retentions. If the TK-Snap or TK-Fric system is only to be used prophylactically, the receiving cavity is closed with the TK-Snap inactive. If friction decreases, a TK-Snap or TK-Fric element can be used as an alternative.

10 When using the TK-Fric system, three different friction strengths are available to control the friction:

- the white element with light friction
- the yellow element with medium friction
- the red element with strong friction

TK-Fric components:

- TK-Fric white, light friction
- TK-Fric yellow, medium friction
- TK-Fric red, strong friction
- Secondary part made of plastic, combustible without residue
- Secondary part HSL, cast-on up to 1400° C to all precious metal alloys (not for non-precious alloys)
- Secondary part Co-Cr-Mo for one-piece casting technique
- TK-Snap inactive

Aids:

- Primary part for the fabrication of the receiving recess in the primary crown
- Primary part for the fabrication of resilience telescopes
- Placeholder for the one-piece casting technique

Tools:

- Parallel holder
- Centering device
- TKR milling cutter for finishing the receiving recess

TK-Fric sets:

- Starter set secondary part made of plastic, combustible without residue
- Starter set secondary part HSL, cast-on up to 1400° C to all precious metal alloys
- Starter set secondary part Co-Cr-Mo for the one-piece casting technique

2. Intended purpose

The TK-Fric retention element manufactured by Si-tec is used with telescopic crowns and other parallel milled dental structures to secure a defined retention force between primary parts attached in the mouth and a removable secondary denture.

3. Indication

Use as a retention element for parallel milled telescopic crowns, bars and custom structures. These can be both periodically and implant-supported.

4. Contraindication

TK-Fric is not suitable for:

- conically milled primary parts.
- severely periodontally damaged abutment teeth.
- patients who, for health reasons, are unable to keep the regularly required monitoring appointments.

TK-Fric must not be used in patients with allergies or suspected allergies to the materials contained in the product.

5. Expected clinical benefit

Protection of the periodontium through defined pull-off forces with stable and comfortable fixation of the denture, which improves the oral health and quality of life of the patient. The clicking when inserting signals the secure fit of the denture to the patient.

6. Target patient group

Patients with residual dentition and completed bone growth for whom fully fixed dentures are not possible or indicated.

7. Intended users

Healthcare professionals (dentists, dental technicians).

The dentist is responsible for explaining to the patients the correct handling during insertion and removal of the denture as well as the correct cleaning.

8. Qualification of healthcare professionals

It is assumed that professional knowledge in professional dentistry or dental technology is available. The current instructions for use must be available and observed. The production and regular check of dentures may only be carried out by qualified professionals.

9. Storage

Keep away from sunlight, in the original packaging.

10. Materials

Texin, HSL (Au, Pt, Pd, Ir), CoCrMo,

11. Side effects / Warnings

Side effects are excluded if used properly.

12. Single use

The product is intended for single-use. The product cannot be re-used because it is firmly connected to the denture and removal leads to a lack of function or destruction of the products.

13. Joint use with other products

For secure fixation of the product in prosthetics, only CE-certified attachment composites (adhesives) approved for the dental sector may be used. According to the manufacturer's instructions, the adhesive used must be explicitly indicated for bonding with the present materials (see point 10). The respective processing and safety instructions of the adhesive manufacturer must be strictly observed.

14. Sterility status

The products are delivered non-sterile.

15. Cleaning / Disinfection

VAH-listed, material-compatible cleaning agents and disinfectants may be used for cleaning and disinfecting the medical device. The following processes are suitable:

- Immersion bath with suitable cleaning agents and disinfectants (e.g. based on quaternary ammonium compounds, amines or peroxides)
- alcoholic immersion bath and spray disinfection ($\approx 70\text{--}80\%$ ethanol/propanol)

The concentrations and exposure times specified by the respective manufacturer must always be observed.

Residues must be removed after cleaning and disinfection with running water of drinking water quality.

16. Reporting of incidents

Serious incidents in connection with the product must be reported by the user to the manufacturer and to the competent authority of the Member State in accordance with the applicable legal requirements.

17. Symbols and their meaning

	Manufacturer		Consult instructions for use
	Date of manufacture		Medical Device
	Batch code		Unique Device Identifier
	Catalogue number		Health Industry Bar Code Standard
	Do not re-use		CE marking
	Keep away from sunlight		

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 info@si-tec.de

 www.si-tec.de

Certified Quality Management System
EN ISO 13485

CE 0044

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**Si-tec GmbH**