TK-Snap rund

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
TK-Snap round

Your Advantages at a Glance:

- secure holding force
- Audible click when inserting
- Secure fixed position even with short friction surfaces
- Data set available for digital processing (additive methods)


Processing Instructions

TK-Snap round


►1

For the TK-Snap system, the primary part forming the retention-groove is put into position using the parallel-holder and the wax-up of the primary crown is completed. The shape from occlusal to approximal surface must be rounded to allow easy sliding of the TK-Snap element.


►2

The primary telescopic crown can now be invested and cast. The retention pins of the primary part used for the parallel-holder can remain on the crown for casting; it is easier to remove them in the metal.

►3

The engagement depth of the TK-Snap round element is 0.3 mm. If this dimension is reduced due to reworking, the bore depth can be restored using the TKR bur.

►4

Lift-off technique


To match the retention groove for the snap-ball, the secondary part is picked with the centering tool and enclosed into the second-part telescope design.

The hole of secondary part is covered with a wax sheet. This prevents liquid wax from flowing uncontrollably into the secondary part.

►5


After casting, all investment material residues must be removed from the secondary part to ensure the secure fit of the TK-Snap.

►6


One-piece casting technique

When using the one-piece casting technique, the placeholder is fixed with wax onto the primary crown so that the centering button fits exactly in the bore of the primary crown. A duplicate model is then created, on which the secondary crown and cast framework is created. After casting, the secondary part is glued into the recess using Si-tec adhesive.

►7

After completion of the denture, the TK-Snap round functional element finally is inserted into the secondary part.

TK-Snap round

Ord.Nr.

Parallel Holder	0122
Centering device TK-Snap round	0206
TKR milling cutter for the post-machining of the receiving recess	0123
Primary part for the preparation of the receiving recess in the primary crown	0101
Secondary part round made of residue-free combustible plastic	0202
Secondary part round HSL, ready for casting up to 1400°C to all precious metal alloys (not for NE alloys)	0242
Secondary part round titanium, for the one-piece-casting technique	0252
Placeholder for one-piece-casting technique	0253
TK-Snap round with titanium ball	0203
TK-Snap round Full plastic	0204
TK-Snap round inactive	0205
Si-tec Adhesive for bonding of the secondary parts in one-piece-casting technique	0850

TK-Snap round

1. Intended purpose

The TK-Snap round retention element manufactured by Si-tec is used with telescopic crowns and other parallel milled dental structures to secure a defined retention force between primary parts attached in the mouth and a removable secondary denture.

2. Indication

Use as a retention element for parallel milled telescopic crowns, bars and custom structures. These can be both periodically and implant-supported.

3. Contraindication

TK-Snap round is not suitable for:

- conically milled primary parts.
- severely periodontally damaged abutment teeth.
- patients who, for health reasons, are unable to keep the regularly required monitoring appointments.

TK-Snap round must not be used in patients with allergies or suspected allergies to the materials contained in the product.

4. Expected clinical benefit

Protection of the periodontium through defined pull-off forces with stable and comfortable fixation of the denture, which improves the oral health and quality of life of the patient. The clicking when inserting signals the secure fit of the denture to the patient.

5. Target patient group

Patients with residual dentition and completed bone growth for whom fully fixed dentures are not possible or indicated.

6. Intended users

Healthcare professionals (dentists, dental technicians).

The dentist is responsible for explaining to the patients the correct handling during insertion and removal of the denture as well as the correct cleaning.

7. Qualification of healthcare professionals

It is assumed that professional knowledge in professional dentistry or dental technology is available. The current instructions for use must be available and observed. The production and regular check of dentures may only be carried out by qualified professionals.

8. Storage

Keep away from sunlight, in the original packaging.

9. Materials

Hostaform MT8U01, CoCrMo, Titanium Grade 4, HSL (Au, Pt, Pd, Ir)

10. Side effects / Warnings

Side effects are excluded if used properly.

11. Single use

The product is intended for single-use. The product cannot be re-used because it is firmly connected to the denture and removal leads to a lack of function or destruction of the products.

12. Joint use with other products

For secure fixation of the product in prosthetics, only CE-certified attachment composites (adhesives) approved for the dental sector may be used. According to the manufacturer's instructions, the adhesive used must be explicitly indicated for bonding with the present materials (see point 9). The respective processing and safety instructions of the adhesive manufacturer must be strictly observed.

13. Sterility status

The products are delivered non-sterile.

14. Cleaning / Disinfection

VAH-listed, material-compatible cleaning agents and disinfectants may be used for cleaning and disinfecting the medical device. The following processes are suitable:

- Immersion bath with suitable cleaning agents and disinfectants (e.g. based on quaternary ammonium compounds, amines or peroxides)
- alcoholic immersion bath and spray disinfection ($\approx 70-80\%$ ethanol/propanol)

The concentrations and exposure times specified by the respective manufacturer must always be observed.

Residues must be removed after cleaning and disinfection with running water of drinking water quality.

15. Reporting of incidents

Serious incidents in connection with the product must be reported by the user to the manufacturer and to the competent authority of the Member State in accordance with the applicable legal requirements.

16. Symbols and their meaning

	Manufacturer		Consult instructions for use
	Date of manufacture		Medical Device
	Batch code		Unique Device Identifier
	Catalogue number		Health Industry Bar Code Standard
	Do not re-use		CE marking
	Keep away from sunlight		

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info@si-tec.de

www.si-tec.de

Certified Quality Management System
EN ISO 13485

CE 0044

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Si-tec GmbH