TK-Soft Ceram

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
CE 0044

TK-Soft Ceram

Precision-adjustable friction element with zirconia screw for all metal-free, parallel-walled dental constructions

TK-Soft Ceram

Your Advantages at a Glance:

- MCS-safe, metal- and allergen-free
- Ideal for metal-free removable dentures
- STL and DME datasets for digital processing available free of charge
- Self-retaining via dowel effect without bonding
- Permanent activation via integrated zirconia screw

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

TK-Soft Ceram


►1

If the digital placeholder for TK-Soft is not yet available in the library, it can be requested as an STL or DME file free of charge. The placeholder is extended toward the primary crown and basally to allow for gap-free positioning. The visible light reflection lines serve for positioning and indicate the original size of the TK-Soft element.


►2

The primary and secondary structures are constructed in the usual manner. A bar stub or pontic connector of appropriate size can be created to accommodate the TK-Soft element.

►3

Select the virtual placeholder from the library and click on "Remove" and "Allow all changes" in the program menu. The digital placeholder is positioned parallel to the telescope axis on the primary part. Its occlusal reflection line must tangentially touch the primary part. The basal reflection line must fit the parallel area above the chamfer.

►4

As the digital placeholder is correctly positioned, it can be punched out using the "Apply" button. Make sure to activate the button "allow any changes", so that at the secondary crown the gap for the friction lamella is punched out. The same procedure applies when using additional friction elements. After the CAD design is completed, it can be passed on to CAM production.

►5

The finished milled, metal-free secondary structure includes the receiving channel for the TK-Soft element.

►6

The fit of the receiving channel can best be checked using the placeholder from the analog technique, as it has exactly the same dimensions as the active element. If the receiving channel is too tight after CAM production, the "+2%" larger digital placeholder can be used for future constructions.

►7

To insert the TK-Soft element into the receiving channel, it is inserted into the mounting aid with the activation screw pointing downward. In this way, it can be safely slid into the channel in the secondary structure. It is not bonded but held securely in the channel via a self-retaining dowel effect.


►8

If necessary, the friction of the telescopic crown can be adjusted by slightly turning the activation screw. The activation path of the lamella is 0.45 mm.

TK-Soft Ceram

	Ord.Nr.		Ord.Nr.
TK-Soft Ceram Starter-Set	2400	Individual parts	
6x Friction elements		Friction element	2401
6x Placeholders		Placeholder	2005
1x Placeholder remover		Placeholder Remover	2022
1x Activating instrument		Mounting aid TK-Soft	2023
1x Mounting aid		Activating instrument	2024

TK-Soft Ceram

1. Intended purpose

The TK-Soft Ceram friction element manufactured by Si-tec serves, in the case of telescopic crowns and other parallel-milled, metal-free dental constructions, to secure a defined holding force between primary parts attached in the mouth and a removable secondary dental prosthesis. The friction can be individually adjusted by turning the activating screw and can also be readjusted later.

2. Indication

Use as a friction element of parallel-milled metal-free telescopic crowns, bars and individual constructions. These can be both periodontally and implant-supported.

3. Contraindication

TK-Soft Ceram is not suitable for:

- conically milled primary parts.
- severely periodontally damaged abutment teeth.
- patients who, for health reasons, are unable to keep the regularly required monitoring appointments.

TK-Soft Ceram must not be used in patients with allergies or suspected allergies to the materials contained in the product.

4. Expected clinical benefit

Protection of the periodontium through defined pull-off forces with stable and comfortable fixation of the denture, which improves the oral health and quality of life of the patient. If necessary, the withdrawal force can be readjusted at any time.

5. Target patient group

Patients with residual dentition and completed bone growth for whom fully fixed dentures are not possible or indicated.

6. Intended users

Healthcare professionals (dentists, dental technicians).

The dentist is responsible for explaining to the patients the correct handling during insertion and removal of the denture as well as the correct cleaning.

7. Qualification of healthcare professionals

It is assumed that professional knowledge in professional dentistry or dental technology is available. The current instructions for use must be available and observed. The production and regular check of dentures may only be carried out by qualified professionals.

8. Storage

Keep away from sunlight, in the original packaging.

9. Werkstoffe

Hostaform MT8U01, zirconium oxide (ZrO₂)

10. Side effects / Warnings

Side effects are excluded if used properly.

11. Single use

The product is intended for single-use. The product cannot be re-used because it is firmly connected to the denture and removal leads to a lack of function or destruction of the products.

12. Joint use with other products

For secure fixation of the product in prosthetics, only CE-certified attachment composites (adhesives) approved for the dental sector may be used. According to the manufacturer's instructions, the adhesive used must be explicitly indicated for bonding with the present materials (see point 9). The respective processing and safety instructions of the adhesive manufacturer must be strictly observed.

13. Sterility status

The products are delivered non-sterile.

14. Cleaning / Disinfection

VAH-listed, material-compatible cleaning agents and disinfectants may be used for cleaning and disinfecting the medical device. The following processes are suitable:

- Immersion bath with suitable cleaning agents and disinfectants (e.g. based on quaternary ammonium compounds, amines or peroxides)
- alcoholic immersion bath and spray disinfection (≈ 70–80 % ethanol/propanol)

The concentrations and exposure times specified by the respective manufacturer must always be observed.

Residues must be removed after cleaning and disinfection with running water of drinking water quality.

15. Reporting of incidents

Serious incidents in connection with the product must be reported by the user to the manufacturer and to the competent authority of the Member State in accordance with the applicable legal requirements.

16. Symbols and their meaning

	Manufacturer		Consult instructions for use
	Date of manufacture		Medical Device
	Batch code		Unique Device Identifier
	Catalogue number		Health Industry Bar Code Standard
	Do not re-use		CE marking
	Keep away from sunlight		

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Certified Quality Management System
EN ISO 13485

CE 0044

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Si-tec GmbH