
Quick-tec System

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

CE 0044

Quick-tec System


For restoring perfect friction in telescopic crowns

Blue -
normal friction

Green -
increased friction

Gold -
normal friction

Your Advantages at a Glance:

- Two friction levels to choose from
- Simple and quick working procedure
- Secure fixation using Si-tec adhesive
- Easy replacement possible
- Low-cost alternative to manufacturing new telescopic crowns

Processing Instructions

Quick-tec System

►1

The required milling in the telescopic crown is done using the two included grinding tools. It is recommended to start with the diamond grinder and then finish the sharp-edged cavity with the carbide cutter. The correct depth of the recess is reached when the rotating shaft produces a shiny spot on the metal. Ensure sufficient cooling during milling.

►2

Since the retention recess is only 0.9 mm deep in the center and tapers off toward the edge, the connection between secondary crown and denture base is not weakened. Milling through the crown wall into the resin area poses no problem, as the recess is closed again by the Quick-tec element.

Now select according to the desired friction:
Quick-tec blue/gold for normal friction or Quick-tec green for increased friction.

►3

Apply adhesive to the dry and grease-free retention recess. Use the adapter to securely pick up the friction element and place it into the recess. Remove the adapter cervically and align the Quick-tec into its correct position.

During the curing process, the Si-tec adhesive forms a watery oxygen inhibition layer on the surface. This can be wiped off with a paper towel after curing.

Remove any remaining adhesive residue manually.


►4

Before the first insertion of the denture, coat the Quick-tec elements with Vaseline. This allows the elements to glide more easily along the primary crown. If the resulting friction is too strong, the central pin of the Quick-tec element can be carefully reduced using a rubber polisher.

Quick-tec System

Ord.Nr.

Ord.Nr.

	Quick-tec Starter-Set 4x Friction elements blue, normal friction 2x Friction elements green, increased friction 1x Adapter 1x Diamant grinder 1x Hard metal milling cutter	0600T	Quick-tec Starter-Set gold 6x Friction elements gold, normal friction 1x Adapter 1x Diamant grinder 1x Hard metal milling cutter	0670T
	Individual parts Friction elements blue, normal friction Friction elements green, increased friction Friction elements gold, normal friction Adapter Diamant grinder Hard metal milling cutter Si-tec Adhesive (for gluing the friction elements)	0601 0602 0603 0621 0623 0624 0850		

Quick-tec system

1. Intended purpose

The Quick-tec manufactured by Si-tec is a repair element for telescopic crown technique.

The Quick-tec repair element is able to restore lost holding forces via friction by protruding into the internal crown space of the secondary telescope. The system is used for subsequent integration into telescopic crowns. It is bonded into a recess to be fabricated with a defined grinding tool. The friction can be individually selected using two different friction elements.

2. Indication

Use as a repair element for parallel-milled telescopic crowns which have insufficient holding force.

3. Contraindication

Quick-tec is not suitable in case of:

- conically milled primary parts.
- patients who, for health reasons, are unable to keep the regularly required follow-up appointments.

In patients with allergies or suspected allergies to the materials contained in the product or the adhesive used, Quick-tec must not be used..

4. Expected clinical benefit

Restoration of the holding force of a telescopic denture, which improves the oral health and quality of life of patients.

5. Target patient group

Patients with residual dentition and completed bone growth for whom fully fixed dentures are not possible or indicated.

6. Intended users

Healthcare professionals (dentists, dental technicians).

The dentist is responsible for explaining to the patients the correct handling during insertion and removal of the denture as well as the correct cleaning.

7. Qualification of healthcare professionals

It is assumed that professional knowledge in professional dentistry or dental technology is available. The current instructions for use must be available and observed. The production and regular check of dentures may only be carried out by qualified professionals.

8. Storage

Keep away from sunlight, in the original packaging.

9. Materials

Texin RxT65D

10. Side effects / Warnings

Side effects are excluded if used properly.

11. Single use

The product is intended for single-use. The product cannot be re-used because it is firmly connected to the denture and removal leads to a lack of function or destruction of the products.

12. Joint use with other products

For secure fixation of the product in prosthetics, only CE-certified attachment composites (adhesives) approved for the dental sector may be used. According to the manufacturer's instructions, the adhesive used must be explicitly indicated for bonding with the present materials (see point 9). The respective processing and safety instructions of the adhesive manufacturer must be strictly observed.

13. Sterility status

The products are delivered non-sterile.

14. Cleaning / Disinfection

VAH-listed, material-compatible cleaning agents and disinfectants may be used for cleaning and disinfecting the medical device. The following processes are suitable:

- Immersion bath with suitable cleaning agents and disinfectants (e.g. based on quaternary ammonium compounds, amines or peroxides)
- alcoholic immersion bath and spray disinfection ($\approx 70-80\%$ ethanol/propanol)

The concentrations and exposure times specified by the respective manufacturer must always be observed.

Residues must be removed after cleaning and disinfection with running water of drinking water quality.

15. Reporting of incidents

Serious incidents in connection with the product must be reported by the user to the manufacturer and to the competent authority of the Member State in accordance with the applicable legal requirements.

16. Symbols and their meaning

	Manufacturer		Consult instructions for use
	Date of manufacture		Medical Device
	Batch code		Unique Device Identifier
	Catalogue number		Health Industry Bar Code Standard
	Do not re-use		CE marking
	Keep away from sunlight		

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info@si-tec.de

www.si-tec.de

Certified Quality Management System
EN ISO 13485

CE 0044

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Si-tec GmbH