
Quick-tec plus

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
CE 0044

Quick-tec plus

The repair element for telescopic crowns,
RS attachments and milling wing


Quick-tec plus

Your Advantages at a Glance:

- Minimal space requirements
- Usable without adhesive
- Simple and quick installation
- Optimal friction behavior
- Long-term stable friction element
- Also suitable for vestibular placement


Processing Instructions

Quick-tec plus


►1

For the friction element positioning in the approximal area is recommended. When using the Quick-tec plus element in an closed dental arch, the ideal position is at an angle of approx. 45° to the oral-vestibular axis. It is important that the Quick-tec plus element is placed towards a parallel surface.


►2

At the selected position, an approach is pre-drilled with the centering drill. This facilitates the secure placement of the carbide twist drill. The telescope wall thickness must not be less than 0.5 mm.


►3

Now, at the marked spot, a receiving hole is drilled into the secondary part using the carbide twist drill. Remaining metal shavings must be removed inside and out with a rubber polisher.


►4

The ball shaped countersink-grinder is then used to grind a trough into the inner wall of the secondary telescope. The correct depth is reached when the shaft touches both sides and leaves a shiny mark.


►5

For wall thicknesses over 0.6 mm, a corresponding reduction must be made from the outside using the countersink grinder to achieve a wall thickness of 0.5 to 0.6 mm.


►6

Insert the Quick-tec plus element from inside to outside through until it clicks into place and rests against the inner wall of the crown.

►7

Cut off the outwarded appendix. A small part of the retention thickening must remain to ensure the Quick-tec plus element is securely held in position.


►8

The result: a barely visible bore sealed with a precision friction element. For discreet placement in veneered areas, we recommend the tooth-coloured Quick-tec plus elements.

Quick-tec plus

Best.Nr.

	Quick-tec plus Starter-Set 4x Friction elements transparent 2x Friction elements tooth-colored 1x Countersink grinder 1x Centering drills 2x HM spiral drill	1400T
	Individual parts Friction element transparent Friction element tooth-colored Countersink grinder Centering drills HM spiral drill	1401 1401ZF 1423 1424 1425

Quick-tec plus

1. Intended purpose

The Quick-tec plus manufactured by Si-tec is a repair element for the telescopic crown technique.

The Quick-tec plus repair element is able to restore lost holding forces via friction by protruding into the internal crown space of the secondary telescope. The system is used for subsequent integration into telescopic crowns. It is inserted into a drilled receiving channel and snaps into place automatically.

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3. Contraindication

Quick-tec plus is not suitable in case of:

- conically milled primary parts
- patients who, for health reasons, are unable to keep the regularly required follow-up appointments.

In patients with allergies or suspected allergies to the materials contained in the product, Quick-tec plus must not be used.

4. Expected clinical benefit

Restoration of the holding force of a telescopic denture, which improves the oral health and quality of life of patients.

5. Target patient group

Patients with residual dentition and completed bone growth for whom fully fixed dentures are not possible or indicated.

6. Intended users

Healthcare professionals (dentists, dental technicians).

The dentist is responsible for explaining to the patients the correct handling during insertion and removal of the denture as well as the correct cleaning.

7. Qualification of healthcare professionals

It is assumed that professional knowledge in professional dentistry or dental technology is available. The current instructions for use must be available and observed. The production and regular check of dentures may only be carried out by qualified professionals.

8. Storage

Keep away from sunlight, in the original packaging.

9. Materials

Ateva G 2803

10. Side effects / Warnings

Side effects are excluded if used properly.

11. Single use

The product is intended for single-use. The product cannot be re-used because it is firmly connected to the denture and removal leads to a lack of function or destruction of the products.

12. Joint use with other products

For secure fixation of the product in prosthetics, only CE-certified attachment composites (adhesives) approved for the dental sector may be used. According to the manufacturer's instructions, the adhesive used must be explicitly indicated for bonding with the present materials (see point 9). The respective processing and safety instructions of the adhesive manufacturer must be strictly observed.

13. Sterility status

The products are delivered non-sterile.

14. Cleaning / Disinfection

For cleaning and disinfection of the medical device, VAH-listed, material-compatible cleaning and disinfecting agents may be used.

The following methods are suitable:

- Immersion bath with suitable cleaning and disinfecting agents (e.g. based on quaternary ammonium compounds, amines or peroxides)

The following methods are conditionally suitable:

- Immersion bath with $\approx 70-80\%$ ethanol/propanol with a short exposure time
- Alcoholic spray disinfection ($\approx 70-80\%$ ethanol/propanol)

The concentrations and exposure times specified by the respective manufacturer must always be observed.

For the Quick-tec plus element, it must be noted that prolonged exposure times in highly concentrated alcohol solution can lead to swelling and material changes.

Residues must be removed after cleaning and disinfection with running water of drinking water quality.

15. Reporting of incidents

Serious incidents in connection with the product must be reported by the user to the manufacturer and to the competent authority of the Member State in accordance with the applicable legal requirements.

16. Symbols and their meaning

	Manufacturer		Consult instructions for use
	Date of manufacture		Medical Device
	Batch code		Unique Device Identifier
	Catalogue number		Health Industry Bar Code Standard
	Do not re-use		CE marking
	Keep away from sunlight		

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info@si-tec.de

www.si-tec.de

Certified Quality Management System
EN ISO 13485

CE 0044

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Si-tec GmbH