Safe-tec II

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
Safe-tec II

1. Processing techniques

Lifting technique with telescopic crowns

1 The primary part is waxed up using the parallel holder parallel to the milling axis of the telescopic crowns. It can be adapted to the course of the alveolar ridge or designed to be free-cleansing. Placeholder "D" can be fixed in the receiving opening of the primary part as an orientation aid.

In this way, the final position of the secondary part can be recognised during waxing up, and positioning errors can be avoided.

2 The approach angle of the primary part should deviate slightly (approx. 15°) from the course of the alveolar ridge. This ensures that in the finished denture, the bolt head is lowered distally, allowing the tongue to glide along with minimal interference and simultaneously enabling easy opening by the patient.

3 The secondary structure can now be fabricated on the finished primary parts.

4 Placeholder "F" is inserted into the receiving opening of the primary part and does not need to be insulated when using modelling resin. The middle section of the placeholder corresponds to the size of the latch secondary part.

5 The design of the secondary structure is carried out in the usual manner. For full lingual veneering, the middle section of the placeholder is encased by approx. 2 mm with modelling resin.


6 The placeholder is pulled out of the wax-up with a rotary movement before lifting the secondary structure.

7 After casting, the secondary structure is fitted. The latch secondary part must be fitted into the receiving channel without tension.

8 For bonding, prepare the surfaces and protect the latch axis with Vaseline against the ingress of adhesive. Now bond the secondary part. The latch axis is seated in the receiving opening of the primary part.

9 During finishing, design the acrylic area around the bolt head such that a mesial recess is available for opening the latch. When closed, the bolt head should not be elevated in the occlusal or distal directions.

Crowns with RS milling and interlock

1 The positioning of the primary part is carried out parallel to the milling axis of the RS milling, as is also the case with telescopic crowns.


2 The finished parallel-milled crowns with milling wing and interlock. In free-end situations, interlock milling is mandatory to compensate for leverage forces.

3 Placeholder "F" is inserted into the receiving opening of the primary part and does not need to be insulated when using modelling resin. The middle section of the placeholder corresponds to the size of the latch secondary part.

4 The design of the secondary structure is carried out in the usual manner. For full lingual veneering, the middle section of the placeholder is encased by approx. 2 mm with modelling resin.

5 The placeholder is pulled out of the wax-up with a rotary movement before lifting the secondary structure.


6 After casting, the secondary structure is fitted. The latch secondary part must be fitted into the receiving channel without tension. The further processing steps are the same as for the fabrication of telescopic crowns.


Replacement of the functional part

Should it be necessary to replace or remove the functional part (inner sleeve with axis), this can be done by unscrewing it using the screw instrument. This is guided around the bolt head and can then engage both slots of the white housing with its blades. The functional part can now be unscrewed. The bonded housing remains in the denture.

Regarding assembly: Insert the functional part into the sleeve and carefully screw it in using the screw instrument. The friction force of the axis can be influenced by the strength of the torque.

Digital fabrication with telescopic crowns (Additive processes)

1 Import the available files for the primary part and the placeholder into the library folders. For subsequent use, please note that the primary part is an additive element and the placeholder is a subtractive element.


2 The positioning of the digital primary part is carried out parallel to the milling axis of the primary telescopes, as in analogue technology.

3 The milled primary parts are scanned and prepared for the secondary structure.

4 The digital placeholder is positioned as a subtractive element in the bore for the latch axis.

5 The design of the secondary structure is carried out in the usual manner. For full lingual veneering, the middle section of the placeholder is encased by approx. 2 mm.

6 The placeholder is removed, and the secondary structure can now be transferred to the CAM area for production. The further processing steps are the same as for analogue processing.

Safe-tec II Components

Main components

Safe-tec II Secondary part, titanium / plastic
(Titanium sleeve, plastic housing, axis)

Safe-tec II Functional part, titanium / plastic
(Plastic housing, axis)

Aids

Primary part, combustible without residue

Placeholder D + F

Tools

Parallel holder

Screwing instrument

Sets

Safe-tec II Initial equipment:

2 Latches complete with aids and tools


Safe-tec II Single pack

1 Latch complete with aids

A suitable adhesive is required for bonding the secondary part (see section 13).

2. Intended purpose

The Safe-tec II bolt manufactured by Si-tec is used for the secure attachment of removable dentures to fixed restorations. These can be designed both unilaterally and bilaterally. The secondary part is fixed in the suprastructure by simple bonding. The required friction force (opening force) of the bolt axis can be adjusted individually according to the needs of the patient.

3. Indication

Use as a connecting element for parallel-milled telescopic crowns, bars and individual constructions with removable dentures. The primary construction can be both periodontally supported and implant-supported.

4. Contraindication

Safe-tec II is not suitable for:

- Patients with impaired fine motor skills.
- Severely periodontally damaged abutment teeth.
- Patients who, for health reasons, are unable to attend the necessary regular follow-up appointments.

Safe-tec II must not be used in patients with allergies or suspected allergies to the materials contained in the product or the adhesive used.

5. Expected clinical benefit

Protection of the periodontium, as no withdrawal forces occur. The denture is securely fixed by the locking mechanism. This serves the oral health and quality of life of the patients. In the case of unilateral use, the dental prosthesis can be reduced to a minimum.

6. Patient target group

Patients with residual dentition and completed bone growth for whom fully fixed restorations are not possible or indicated.

7. Intended users

Healthcare professionals (dentists, dental technicians).

The dentist is responsible for explaining to the patients the correct handling during the insertion and removal of the denture as well as the correct cleaning.

8. Qualification of healthcare professionals

It is assumed that specialist knowledge in professional dentistry or dental technology is available. The current instructions for use must be available and observed. The manufacture and regular inspection of dental prostheses may only be carried out by qualified specialists.

9. Storage

Protected from sunlight, in the original packaging.

10. Materials

Tecaform AH, Titanium Grade 4.

11. Side effects / Warnings

Side effects are excluded if used correctly.

12. Single-use

The product is intended for single-use only. The product cannot be re-used because it is bonded to the denture and removal leads to loss of function or destruction of the products.

13. Use in combination with other products

For the secure fixation of the product in prosthetics, only CE-marked and dental-approved luting composites (adhesives) may be used. According to the manufacturer's instructions, the adhesive used must be explicitly indicated for the bonding with the

present materials (see section 9). The respective processing instructions and safety instructions of the adhesive manufacturer must be strictly observed.

14. Sterile state

The products are delivered non-sterile.

15. Cleaning / Disinfection

VAH-listed, material-compatible cleaning and disinfecting agents may be used for cleaning and disinfecting the medical device. The following procedures are suitable:

- Immersion bath with suitable cleaning and disinfecting agents (e.g. based on quaternary ammonium compounds, amines or peroxides).
- Alcoholic immersion bath and spray disinfection (≈ 70 – 80% ethanol/propanol).

The concentrations and contact times specified by the respective manufacturer must always be observed.

Residues must be removed after cleaning and disinfection with running water of drinking water quality in accordance with the TrinkwV (Drinking Water Ordinance).

16. Reporting of incidents

Serious incidents in connection with the product must be reported by the user in accordance with the applicable legal requirements to the manufacturer and to the competent authority of the Member State.

17. Symbols and their meaning

	Manufacturer		Consult IFU
	Date of manufacture		Medical device
	Batch code		UDI identification
	Catalogue number		Health Industry Bar Code Standard
	Do not reuse		CE marking
	Keep away from Sunlight		

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Certified Quality Management System
EN ISO 13485

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Si-tec GmbH