
SIM-tec Geschiebe

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
SIM-tec

1. Processing techniques

- 1 The primary part made of combustible without residue plastic is waxed up to the crown close to the stump and parallel to the milled milling wing using the parallel holder.
- 2 When modeling the crown, the surface of the primary part must be kept wax-free.
- 3 After casting, the primary part can be shortened down to 2.9 mm. Round the occlusal surface so that the secondary part is easier to fit. The attachment surfaces are carefully polished in order to remove as little material as possible.
- 4 For the production of the cast framework, shorten the yellow secondary part according to the length of the primary part. Prepare the model for duplicating, while blocking out the basal area under the attachment.

- 5 In the finished cast framework, the secondary part is inserted into the receiving opening using the mounting aid. The retention lug of the secondary part must securely click into place. If necessary, the receiving opening can be carefully reworked.
- 6 The secondary part is securely held by the integrated retention lug and is not bonded.
- 7 Should a higher friction be desired, the yellow secondary part can be replaced with a red secondary part.

SIM-tec components:


- SIM-tec primary part
- SIM-tec secondary part yellow, normal friction
- SIM-tec secondary part red, increased friction

Tools:

- parallel holder
- mounting aid

Sets:

- SIM-tec starter set

2. Intended purpose

The SIM-tec attachment manufactured by Si-tec is a plastic T-attachment for connecting removable dentures with fixed dentures. The primary part consists of residue-free burn-out plastic and can thus be cast in any dental alloy together with the crowns. With the easily replaceable secondary parts in two friction strengths, the holding force of the denture can be selected.

3. Indication

Use as a connecting element between fixed and removable dentures.

4. Contraindication

SIM-tec is not suitable in case of:

- severely periodontally damaged abutment teeth.
- patients who, for health reasons, are unable to keep the regularly required follow-up appointments.

In patients with allergies or suspected allergies to the materials used in the product, SIM-tec must not be used.

5. Expected clinical benefit

Due to the conical bevel, easy insertion of the denture by the patient is possible. Tilting is avoided, which protects the periodontium. This improves the oral health and quality of life of patients.

6. Target patient group

Patients with residual dentition and completed bone growth for whom fully fixed dentures are not possible or indicated.

7. Intended users

Healthcare professionals (dentists, dental technicians).

The dentist is responsible for explaining to the patients the correct handling during insertion and removal of the denture as well as the correct cleaning.

8. Qualification of healthcare professionals

It is assumed that professional knowledge in professional dentistry or dental technology is available. The current instructions for use must be available and observed. The production and regular check of dentures may only be carried out by qualified professionals.

9. Storage

Keep away from sunlight, in the original packaging.

10. Materials

Hostaform MT8U01

11. Side effects / Warnings

Side effects are excluded if used properly.

12. Single use

The product is intended for single-use. The product cannot be re-used because it is firmly connected to the denture and removal leads to a lack of function or destruction of the products.

13. Joint use with other products

For secure fixation of the product in prosthetics, only CE-certified attachment composites (adhesives) approved for the dental sector may be used. According to the manufacturer's instructions, the adhesive used must be explicitly indicated for bonding with the present materials (see point 10). The respective processing and safety instructions of the adhesive manufacturer must be strictly observed.

14. Sterility status

The products are delivered non-sterile.

15. Cleaning / Disinfection

VAH-listed, material-compatible cleaning agents and disinfectants may be used for cleaning and disinfecting the medical device. The following processes are suitable:

- Immersion bath with suitable cleaning agents and disinfectants (e.g. based on quaternary ammonium compounds, amines or peroxides)
- alcoholic immersion bath and spray disinfection ($\approx 70-80\%$ ethanol/propanol)

The concentrations and exposure times specified by the respective manufacturer must always be observed.

Residues must be removed after cleaning and disinfection with running water of drinking water quality.

16. Reporting of incidents

Serious incidents in connection with the product must be reported by the user to the manufacturer and to the competent authority of the Member State in accordance with the applicable legal requirements.

17. Symbols and their meaning

	Manufacturer		Consult instructions for use
	Date of manufacture		Medical Device
	Batch code		Unique Device Identifier
	Catalogue number		Health Industry Bar Code Standard
	Do not re-use		CE marking
	Keep away from sunlight		

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info@si-tec.de

www.si-tec.de

Certified Quality Management System
EN ISO 13485

CE 0044

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Si-tec GmbH