Grazilklammern

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Denture clasps

1. General information

- Simple standard solution for the patient
- Universally applicable
- Manufactured from orthodontic steel
- Fast and uncomplicated processing

2. Intended purpose

The slender clasps manufactured by Si-tec are prefabricated retaining clasps suitable for the manufacture of removable dentures. They can be individually adapted to the tooth shape if required.

► denture clasp

3. Indication

Use as a simple retaining element for dentures on remaining dentition. This technique enables particularly easy insertion and removal of the denture.

4. Contraindication

In patients with allergies or suspected allergies to the materials contained in the product, the slender clasps must not be used.

5. Expected clinical benefit

Secure retention of the denture, which improves the oral health and quality of life of patients. The denture can be manufactured without pre-treatment or processing of the remaining dentition.

6. Target patient group

Patients with residual dentition and completed bone growth for whom fully fixed dentures are not possible or indicated.

► denture clasp

7. Intended users

Healthcare professionals (dentists, dental technicians).

The dentist is responsible for explaining to the patients the correct handling during insertion and removal of the denture as well as the correct cleaning.

8. Qualification of healthcare professionals

It is assumed that professional knowledge in professional dentistry or dental technology is available. The current instructions for use must be available and observed. The production and regular check of dentures may only be carried out by qualified professionals.

9. Storage

Keep away from sunlight, in the original packaging.

10. Material

Stainless steel 1.4310

11. Side effects / Warnings

Side effects are excluded if used properly.

12. Single use

The product is intended for single-use. The product cannot be re-used because it is firmly connected to the denture and removal leads to a lack of function or destruction of the products.

13. Joint use with other products

For secure fixation of the product in prosthetics, only CE-certified attachment composites (adhesives) approved for the dental sector may be used. According to the manufacturer's instructions, the adhesive used must be explicitly indicated for bonding with the present materials (see point 10). The respective processing and safety instructions of the adhesive manufacturer must be strictly observed.

14. Sterility status

The products are delivered non-sterile.

15. Cleaning / Disinfection

VAH-listed, material-compatible cleaning agents and disinfectants may be used for cleaning and disinfecting the medical device. The following processes are suitable:

- Immersion bath with suitable cleaning agents and disinfectants (e.g. based on quaternary ammonium compounds, amines or peroxides)
- alcoholic immersion bath and spray disinfection ($\approx 70\text{--}80\%$ ethanol/propanol)

The concentrations and exposure times specified by the respective manufacturer must always be observed.

Residues must be removed after cleaning and disinfection with running water of drinking water quality.

16. Reporting of incidents

Serious incidents in connection with the product must be reported by the user to the manufacturer and to the competent authority of the Member State in accordance with the applicable legal requirements.

17. Symbols and their meaning

	Manufacturer		Consult instructions for use
	Date of manufacture		Medical Device
	Batch code		Unique Device Identifier
	Catalogue number		Health Industry Bar Code Standard
	Do not re-use		CE marking
	Keep away from sunlight		